



Dentalone

The dentist's office Mobile cart



Army & Hospitals



Humanitarian



Home Care



Schools





NSK Dentabone

CONTROL PANEL

UP

DOWN

Spray Volume Scal

Spray Volume NLX

Spray ON OFF

UAL LOCK

LEVEL

EXTENDED FUNCTIONALITIES

ADVANTAGES

- Compact
- Simple to use
- Mobile
- Plug and play system
- Easy to carry
- User-friendly
- Extensible to Cart function
- Large range functionalities

Dentalone is the dentist's office mobile cart easily expandable from "carry" position to a comfortable, user-friendly dental treatment cart ready to use anywhere, anytime.

Dentalone offers a large range of instruments and extended functionalities, thus practitioners can realise dental treatments outside any dentist's office:

- Instruments carrier with automatic selection of instruments
- Illuminated LCD control panel
- Brushless NLX nano micromotor with LED Light
- Integrated Multi Function Ultrasonic Scaler Varios 170 with LED Light
- Syringe with air/water/spray with LED Light
- Cannula for saliva

DENTIST'S OFFICE MOBILE CART

Plug and play



User-friendly



NSK

Dentalone

200000 ^{min⁻¹}
NLX nano
PRG HIGH GEAR 1:5
1

FWD/REV

MODE

SELECT

SET UP

ESC

ENT/SAVE

NSK

Varios2 Lux

NSK

Ti-Max Z95L

NSK

NLX nano

NSK

FEATURES

CONTENT

- Carrying case
- Water tank for spray
- tank for vaccum
- Built in compressor
- Built in suction pump
- Variable foot-control
- Extendable tablet

SPECIFICATIONS

Dimensions (HxWxD)

- | | |
|--|----------------------|
| • Case with trolley (unpulled trolley handle): | 545 x 351 x 405 mm |
| • Case with trolley (pulled trolley handle): | 1,075 x 351 x 405 mm |
| • Dimensions without case cover: | 496 x 350 x 405 mm |
| • Working height: | 766 mm |
| • Weight with case: | 18,5 kg |
| • Power Supply: | 100–230 V~ 50/60 Hz |

Compact



Expandable to Cart Function

496 mm



766 mm





EC REP

NAKANISHI INC.

www.nsk-inc.com

700 Shimohinata
Kanuma-shi
Tochigi 322-8666,
Japan

TEL : +81 (0)289-64-3380
FAX : +81 (0)289-62-5636

NSK Europe GmbH

www.nsk-europe.de

Elly-Beinhorn-Strasse 8
65760 Eschborn,
Germany

TEL : +49 (0)6196 77606-0
FAX : +49 (0)6196 77606-29

NSK United Kingdom Ltd

www.nsk-uk.com

Office 5, Gateway1000,
Arlington Business Park,
Whittle Way,
Stevenage, SG1 2FP, UK
TEL : +44 (0)1438 310670
FAX : +44 (0)1438 310671

NSK Dental LLC

www.nskdental.com

1800 Global Parkway
Hoffman Estates, IL 60192,
USA

TEL : +1 888 675 1675
FAX : +1 800 838 9328

NSK Oceania Pty Ltd

www.nskoceania.com.au

Unit 22, 198-222 Young St.
Waterloo, Sydney,
NSW 2017, Australia

TEL : +61 (0)2 8306 3000
FAX : +61 (0)2 8306 3030

NSK Middle East

www.nsk-inc.com

Room 6EA-701, 7th Floor,
East Wing No.6
Dubai Airport Free Zone,
PO Box 54316 Dubai, UAE
TEL : +971 (0)4-204-5507
FAX : +971 (0)4-204-5508

NSK Asia Pte Ltd

www.nsk-inc.com

1 Maritime Square,
#09-33 HarbourFront Centre,
Singapore 099253

TEL : +65-6376-5018
FAX : +65-6376-0233

NSK France

www.nsk.fr

32, rue de Lisbonne
75008 Paris
France

TEL : +33(0)1 567 959 80
FAX : +33(0)1 567 959 81

NSK Dental Spain SA

www.nsk-spain.es

C/ M6dena, 43 El Soho - Eur6polis
28232 Las Rozas, Madrid - Espa1a
TEL : +34-91-626-61-28
FAX : +34-91-626-61-32

Specifications are subject to change without notice.