




CREATE IT.


TIP




AIR SCALER TIPS & APPLICATIONS

GUIDE


A Tip range to suit your needs





SCALING



SCALING



SCALING

MODEL S1

- Removal of Supra gingival calculus.
- Rimval of marginal gingival calculus.





ORDER CODE **Z252411**

• Universal Type

MODEL S2

- Removal of Supra gingival calculus.




ORDER CODE **Z252413**

• Sickle Type



MODEL S3


- Removal of Supra gingival calculus.
- Rimval of marginal gingival calculus.



ORDER CODE **Z252412**

• Perio Type



PERIO



PERIO

MODEL S10

- Removal of sub gingival calculus.
- Removal of narrow interdental space calculus.





ORDER CODE **Z252420**

- Universal Type Slim Line

MODEL S20

- Removal of sub gingival calculus.
- Removal of sub gingival biofilm





ORDER CODE **Z252441**

- Root Planning (Straight tip)

MODEL S21R

- Removal of sub gingival calculus.
- Removal of narrow interdental space calculus.
- Removal of sub gingival biofilm





ORDER CODE **Z252442**

- Root Planning (Right curved tip)

MODEL S21L

- Removal of sub gingival calculus.
- Removal of narrow interdental space calculus.
- Removal of sub gingival biofilm




ORDER CODE **Z252443**

- Root Planning (Left curved tip)




PERIO (Maintenance Tip)

For the maintenance of prostheses, implant-supported prostheses. These tips are made from a special material which will not damage prostheses and restorative materials.



MODEL V-P10

V-P10 is for the maintenance of prostheses, its shape offers easy access to all treatment areas.




ORDER CODE **Y900184**

• Pack of 3 • S35 holder is not included



MODEL V-P11R

V-P11R/L left and right are double bent tips. This tip is extra long to easily access difficult to reach areas and to avoid interference with neighbouring teeth or contours. Its shape was designed to best fit different tooth morphologies, making this a versatile tip for a wide range of uses from anterior teeth to premolars.



ORDER CODE **Y1002165**

• Right curved type • Pack of 3 • S35 holder is not included



MODEL V-P11L




ORDER CODE **Y1002166**

• Left curved type • Pack of 3 • S35 holder is not included

MODEL V-P12

V-P12 is a longer version of the popular V-P10 tip. It serves a wide range of applications.



ORDER CODE **Y1002167**

• Pack of 3 • S35 holder is not included




MODEL S35



ORDER CODE **Z252015**

• Includes E-Tip replacement wrench • Plastic Tip is not included

MODEL Perio-Control Tips Air Scaler Kit




ORDER CODE **Y1002754**

• S35 holder : 1pc.
• Includes E-Tip replacement wrench
• V-P10, V-P11R, V-P11L, V-P12: 1 pc. each

RESTORATIVE



RESTORATIVE (For Minimal Intervention)








MODEL **S65D**

MODEL **S66D**

MODEL **S67D**

MODEL **S68D**

The spherical shape of these tips make them suitable for delicate adjustments as well as caries removal in the cervical and interproximal region.



ORDER CODE **Z252465**

ORDER CODE **Z252466**

ORDER CODE **Z252467**


ORDER CODE **Z252468**

• Diamond coated

• Diamond coated


• Diamond coated

• Diamond coated




Caries removal in the cervical and interproximal region.

RESTORATIVE





RESTORATIVE (For Preparation of crown cores)



MODEL **S71D**

MODEL **S72D**

These tips are safe sided and flat for minimal mesial or distal preparations.



ORDER CODE **Z252471**


ORDER CODE **Z252472**

• Diamond coated

• Diamond coated



RESTORATIVE (For Preparation of crown cores)









MODEL **S81D**

MODEL **S82D**

MODEL **S83D**

MODEL **S86D**

The tips can be used for finishing margins after prepping the area with a conventional dental bur.
That enables you to finish minute sharp edges and gingival sulcus without damage, something that cannot be accomplished with rotary instruments alone.



ORDER CODE **Z252469**

ORDER CODE **Z252462**

ORDER CODE **Z252461**



ORDER CODE **Z252460**


• Diamond coated

• Diamond coated

• Diamond coated

• Diamond coated







V-TIP SYSTEM

V-TIP SYSTEM

NSK V-Tip system is innovative, multi-functional and economical. The V-Tip holder is a separate component to the actual cutting tips and can be easily and inexpensively replaced. NSK V-Tip holder system offers the flexibility and expandability to meet the wide requirements for multi-purpose clinical applications.



MODEL **S75**



ORDER CODE **Z252035**


• Included E-Tip replacement wrench

V-TIP SYSTEM (Endodontics)



NSK offers Endodontic files

MODEL **V-U15**



ORDER CODE **Y900177**

• V-file Irrigation 23 mm #15 • ISO15 • Pack of 6


MODEL **V-U20**



ORDER CODE **Y900178**

• V-file Irrigation 23 mm #20 • ISO20 • Pack of 6


MODEL **V-U25**



ORDER CODE **Y900179**

• V-file Irrigation 23 mm #25 • ISO25 • Pack of 6


MODEL **V-U30**



ORDER CODE **Y900180**

• V-file Irrigation 23 mm #30 • ISO30 • Pack of 6

MODEL **V-U35**



ORDER CODE **Y900181**

• V-file Irrigation 23 mm #35 • ISO35 • Pack of 6



V-TIP SYSTEM

